

# PepBay: Implementation of Bayesian inference in the analysis of peptide arrays


Katarzyna Sidorczuk

University of Wrocław, Faculty of Biotechnology, Department of Bioinformatics  
and Genomics

Why R? 2019, Warsaw

# What are peptide arrays?

- Collections of short protein fragments;
- Efficient tool for search of new biomarkers;
- Peptide array data:
  - very small sample size (patients),
  - large number of variables (peptides),
  - correlated.



# Why traditional methods fail when $p \gg n$ ?

Multiple simultaneous statistical tests  
(control of the 1st type error rate)



Correction for multiple testing




High corrected p-values



Cannot distinguish between noise and significant results

# Solution: Bayesian inference



(Doll, J. C. and Jacquemin, S. J. 2018)

Bayes' theorem:

$$P(\theta|X) = \frac{P(X|\theta) \times P\theta}{P(X)}$$

# Implementation using package BEST

```
y1 <- rnorm(100)  
y2 <- rnorm(100)  
test <- BESTmcmc(y1,y2)
```

```
## Waiting for parallel processing to complete...done.
```

BEST package:

- Based on JAGS;
- Core function:  
BESTmcmc;
- Convenient wrapper:  
tidybayes;
- Alternative: rSTAN.

# Implementation using package BEST

```
y1 <- rnorm(100)
y2 <- rnorm(100)
test <- BESTmcmc(y1,y2)
```

```
## Waiting for parallel processing to complete...done.
```

```
test
```

```
## MCMC fit results for BEST analysis:
## 100002 simulations saved.
##
```

	mean	sd	median	HDIlo	HDIup	Rhat	n.eff
## mu1	-0.04564	0.10282	-0.04571	-0.24666	0.1572	1	58812
## mu2	0.11573	0.09411	0.11573	-0.07092	0.2998	1	60291
## nu	41.47526	30.51176	33.05015	4.66915	102.8904	1	20892
## sigma1	0.98335	0.07812	0.98014	0.83142	1.1381	1	46892
## sigma2	0.90946	0.07362	0.90715	0.76636	1.0552	1	44496

```
##
## 'HDIlo' and 'HDIup' are the limits of a 95% HDI credible interval.
## 'Rhat' is the potential scale reduction factor (at convergence, Rhat=1).
## 'n.eff' is a crude measure of effective sample size.
```

BEST package:

- Based on JAGS;
- Core function: BESTmcmc;
- Convenient wrapper: tidybayes;
- Alternative: rSTAN.

# Implementation using package BEST

BEST package:

- Based on JAGS;
- Core function: BESTmcmc;
- Convenient wrapper: tidybayes;
- Alternative: rSTAN.

```
y1 <- rnorm(100)
y2 <- rnorm(100)
test <- BESTmcmc(y1,y2)
```

```
## Waiting for parallel processing to complete...done.
```

```
test
```

```
## MCMC fit results for BEST analysis:
## 100002 simulations saved.
##          mean      sd    median    HDILO    HDIup Rhat n.eff
## mu1      -0.04564  0.10282 -0.04571 -0.24666  0.1572  1 58812
## mu2       0.11573  0.09411  0.11573 -0.07092  0.2998  1 60291
## nu       41.47526 30.51176 33.05015  4.66915 102.8904  1 20892
## sigma1    0.98335  0.07812  0.98014  0.83142  1.1381  1 46892
## sigma2    0.90946  0.07362  0.90715  0.76636  1.0552  1 44496
##
## 'HDILO' and 'HDIup' are the limits of a 95% HDI credible interval.
## 'Rhat' is the potential scale reduction factor (at convergence, Rhat=1).
## 'n.eff' is a crude measure of effective sample size.
```


```
summary(test)
```

##	mean	median	mode	HDI%	HDILO	HDIup	compVal	%>compVal
## mu1	-0.0456	-0.0457	-0.0470	95	-0.2467	0.157		
## mu2	0.1157	0.1157	0.1089	95	-0.0709	0.300		
## muDiff	-0.1614	-0.1620	-0.1703	95	-0.4375	0.109	0	12.4
## sigma1	0.9833	0.9801	0.9783	95	0.8314	1.138		
## sigma2	0.9095	0.9072	0.9006	95	0.7664	1.055		
## sigmaDiff	0.0739	0.0728	0.0724	95	-0.1213	0.280	0	76.9
## nu	41.4753	33.0502	19.2860	95	4.6691	102.890		
## log10nu	1.5122	1.5192	1.5352	95	0.9187	2.098		
## effSz	-0.1714	-0.1713	-0.1800	95	-0.4670	0.113	0	12.4

# Bayesian inference vs. frequentist methods

Advantages of Bayesian inference:

- complete distributions of reliable values;
- effect size instead of p-value.





# PepBay app screenshots

## Peptide browser

Overview Detailed view n-gram panel

- prot\_id: ID of protein.
- Sequence: a sequence of the peptide. The search box supports partial matching.
- coef\_bin: the phenotype associated with a peptide.
- phens: phenotypes against which the fold difference of measurement is larger than 1.5.
- p-value: raw (non-adjusted) p-value.
- marker\_against: against how many phenotypes the fold difference of measurement is larger than 1.5.

Select peptides individually by clicking on the rows or use the checkbox below to select all.

☐ Select all peptides

Copy CSV Excel Print

prot_id	gene_name	Sequence	coef_bin	p_ben	p_mal	p_WT	marker_against
All	All	All	All	All	All	All	All
NP_000033	APOH	PDNGFVNYPKPK	mal			0.0997	1
NP_000081	COL3A1	EKGSPGAQGPPG	mal			0.0050	1
NP_000081	COL3A1	GEKGEQGPFGVA	mal			0.0046	1
NP_000081	COL3A1	PPGMPPGPRGSPG	mal			0.0026	1


# PepBay app screenshots

## Peptide browser

Overview Detailed view n-gram panel

## Peptide plot

The color of squares represents the phenotype associated with a peptide, triangles represent p-value and boxplot distribution of measured points.



# PepBay app screenshots


## Peptide browser

[Overview](#) [Detailed view](#) [n-gram panel](#)

### n-gram analysis of amino acid motifs in peptides

Motifs are automatically selected by QuiPT (significance level: 0.05). Longer n-grams require longer computation time.

Length of n-grams



# Acknowledgments

- Andreas Weinhäusel  
Austrian Institute of Technology, Department of Molecular Diagnostics
- Michał Burdukiewicz  
Warsaw University of Technology, Faculty of Mathematics and Information Science